

Korus order tracking allows you to quickly track orders, tubes and results for each patient.

Features:

- Check to see if you have orders for which no tubes have been received by the lab
- View orders that are flagged with alerts and exceptions, partial tubes received and cancelled tests
- Filter orders by date range, modality, flag type, and order status

Status And Icons Description

Order Statuses:

- **Not Received** - No tubes for the order have been scanned at the lab.
- **In Process** - At least one tube has been scanned at the lab.
- **Final** - All tests are resulted or cancelled.

Order With Flags:

- **Alerts and Exceptions** - The order has a general lab alert or exception flag.
Note: **Microbiology alerts are not included.**
- **Cancelled Tests** - The order has at least one cancelled test.
- **Partial Tubes Received** - At least one tube for the order has been scanned at the lab and at least one tube still has a status of not received.

- View all patients within your specified date range and drill down to the individual patient requisition, which provides a tube summary and order details

Tube Summary:

- Shows all the tubes that were received by the lab and the time at which they were scanned
- Provides comments regarding specimen condition – e.g., clotted, hemolysis

Order Details:

- Shows all results with alert/exception flags and additional comments from the lab

Tube Summary 4 out of 4 received				
Accession	Tube Type	Specimen	Received	Condition
120481DWL	Lavender	Blood	07/05/2017 11:30 AM PT	
120481DWL	Lavender	Plasma	07/05/2017 11:30 AM PT	
130481DWL	SST Gel	Serum	07/05/2017 12:49 PM PT	Slight Hemolysis
230481DWL	Gold Post_8	Plasma	07/05/2017 16:08 PM PT	

Order Details # of Tests: 35 Frequency: QTR Group By: Report Header					
Test	Result	Units	Reference	Comments	
Hematology					
White Blood Cell Count	6.28	1000/mcL	4.80-10.80		
Red Blood Cell Count	3.80	(L) mill/mcL	4.20-5.40		
Hemoglobin	11.5	(L) g/dL	12.0-16.0		