

Advanced Environmental Reporting

Korus provides unmatched environmental reporting tools to help you monitor the safety of dialysis treatment and comply with the latest quality standards.^{1, 2}

Korus environmental reports help you:

- Document the corrective action steps taken to achieve chemical, bacterial endotoxin and colony count levels that fall within acceptable AAMI-defined ranges
- View the trending of chemical contaminants for a particular source/equipment at both the facility and corporate level, and easily identify your outlier results
- Determine chemical contaminant levels for all equipment in your facility within a specified date range
- Quickly identify all sources/equipment that require immediate action due to outlier bacterial endotoxin, colony count and chemical contaminant levels
 - Sources can be viewed at both the corporate and facility levels

Dialysis Ctr		Facility Summary					
Colony Count: Action Level: 50		Allowable Limit: 200					
Endotoxin: Action Level: 1.00		Allowable Limit: 2.00					
Equipment Name	Test Name	Collected	Days from Last Sample	Source	Result	Units	Action Steps Taken
Bicarb Delivery System (Bds) 1, Return 1	Colony Count	09/06/10 BA*	28	Bicarbonate	<2	CFU/mL	_____
Bicarb Mixer 1, Food	Colony Count	09/06/10 BA*		Water	<2	CFU/mL	_____
Dialysis Machine, Machine 18	Colony Count	09/06/10 BA*	84	Dialysate	<2	CFU/mL	_____
Dialysis Machine, Machine 19	Colony Count	09/06/10 BA*		Dialysate	<2	CFU/mL	_____
Dialysis Machine, Machine 35	Colony Count	09/06/10 BA*	84	Dialysate	<2	CFU/mL	_____
H2o Delivery System (Wds) 1, Return Loop 1	Colony Count	09/06/10 BA*	28	Water	<2	CFU/mL	_____
H2o Delivery System (Wds) 1, Ultrafilter Pool 1	Colony Count	09/06/10 BA*	28	Water	<2	CFU/mL	_____

1. Centers for Medicare and Medicaid Services. 2008; Conditions for Coverage for End-Stage Renal Disease Facilities. 42 CFR Part 494.80 & 494.90.

2. ANSI/AAMI RD62:2001 and ANSI/AAMI RD52:2004. Available at www.aami.org/standards/



Spectra Laboratories, Inc.
 525 Sycamore Drive • Milpitas, CA 95035 • 800-433-3773
 8 King Road • Rockleigh, NJ 07647 • 800-522-4662
www.spectra-labs.com

© 2010 Fresenius Medical Care Holdings, Inc. All rights reserved.
 Spectra, the Spectra logo, Korus, and the Korus logo are trademarks of
 Fresenius Medical Care Holdings, Inc. or its affiliated companies. All other
 trademarks are the property of their respective owners.